3M Infection Prevention Solutions

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United Nations
Environment
Programme

Distr.
GENERAL

UNEP/OzL.Pro/WG.1/21/INF/3
20 June 2001

ORIGINAL: ENGLISH

OPEN-ENDED WORKING GROUP OF THE PARTIES TO
THE MONTREAL PROTOCOL ON SUBSTANCES THAT
DEplete THE OZONE LAYER
Twenty-first meeting
Montreal, 24 June 2001

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1. Following the request from the Technology and Economic Assessment Panel, the Ozone Secretariat wishes to provide Parties with important information likely to help avoid potential confusion over the issue of ethylene oxide sterilization and ozone layer protection.
2. Sterilization is critical to public health and ethylene oxide is used as a very effective sterilizing agent. It can be used in sterilization technology in a mixture with ozone-depleting substances (chlorofluorocarbons (CFC) or hydrochlorofluorocarbons (HCFC)).

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4. Ethylene oxide sterilization presents a potential health and safety hazard in some countries, ethylene oxide use in sterilization technology in order to manage these health and safety hazards. Health and safety hazard considerations for ethylene oxide are independent of considerations of ozone-depleting substances used in ethylene oxide sterilization.
5. Parties are reminded that measures taken under the Montreal Protocol to phase out the use of ozone-depleting substances in ethylene oxide sterilization are taken for the purpose of protecting the ozone layer only.

K0110161 250601

For reasons of economy, this document is printed in a limited number. Delegates are kindly requested to bring their copies to meetings and not to request additional copies.

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
*Ref.: 'Employee Safety and Sterilant Gases', D.J. Hilliker,
Infection Control Today, 1998, pp. 37-45.*

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CDC Home
Centers for Disease Control and Prevention
CDC 24/7: Saving Lives. Protecting People.™

NIOSH
All CDC Topics
Choose a topic above

A-Z Index for All CDC Topics

NIOSH Publications and Products

NIOSH > NIOSH Publications & Products > NIOSH-Issued Publications

Recommend | o | Tweet | Share

DHHS (NIOSH) Publication Number 2005-149 September 2007

NIOSH Pocket Guide to Chemical Hazards

The NIOSH Pocket Guide to Chemical Hazards (NPG) is intended as a source of general industrial hygiene information on several hundred chemicals/classes for workers, employers, and occupational health professionals. The NPG does not contain an analysis of all pertinent data, rather it presents key information and data in abbreviated or tabular form for chemicals or substance groupings (e.g. cyanides, fluorides, manganese compounds) that are found in the work environment. The information found in the NPG should help users recognize and control occupational chemical hazards.

NIOSH Pocket Guide to Chemical Hazards (3rd Printing with minor technical revisions) [PDF 6,192KB]

View a current version that is updated of the [NIOSH Pocket Guide](#).

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Centers for Disease Control and Prevention
National Institute for Occupational Safety and Health

3M Medical Markets Laboratory Europe News

Low Temperature Steam Formaldehyde

World Health Organization

Press Release N° 153

Page 1 of 2

INTERNATIONAL AGENCY FOR RESEARCH ON CANCER

PRESS RELEASE
N° 153

WHO

15 June 2004

IARC CLASSIFIES FORMALDEHYDE AS CARCINOGENIC TO HUMANS

"Twenty-six scientists from 10 countries evaluated the available evidence on the carcinogenicity of formaldehyde, a widely used chemical", reports Dr Peter Boyle, Director of the International Agency for Research on Cancer (IARC), part of the World Health Organization. The working group, convened by the *IARC Monographs Programme*, concluded that formaldehyde is *carcinogenic to humans*. Previous evaluations, based on the smaller number of studies available at that time, had concluded that formaldehyde was *probably carcinogenic to humans*, but new information from studies of persons exposed to formaldehyde has increased the overall weight of the evidence.


The 3M logo, consisting of the letters "3M" in a bold, red, sans-serif font.

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Ref.: 'Mortality Among Workers Exposed TO EO', K. Steenland et al, England Journal of Medicine, May, 1991, pp. 1402-1407.

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Annual WFHSS and JSMI Conference 2012
13th World Sterilization Congress

Problems on Hydrogen Peroxide Sterilisation — New Proposal for Safety and Effective Use —

Rika Yoshida, Hiroyoshi Kobayashi

Division of Infection Prevention and Control,
Tokyo Healthcare University Postgraduate School

23 Nov, 2012

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Background

- Some central service employees of a tertiary university hospital engaged in hydrogen peroxide (HP) sterilisation had experienced the irritations to eyes and throats, and one of them had acute flare and swelling around eyes.
- It was supposed as the adverse event of HP sterilisation as reported to Food and Drug Administration (FDA) in the United States, though all of them were not evaluated as the actual adverse event by HP.

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No regulation for the environmental concentrations
of hydrogen peroxide in Japan

A safe argument has not been held enough
in Japan

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Method 1

- The concentrations of hydrogen peroxide vapour (HPV) in the proximity of HP sterilisers (low-temperature hydrogen peroxide gas plasma steriliser: STERRAD[®] NX[™], STRRAD[®] 100S, and STERRAD[®] 200, Johnson & Johnson, and low-temperature vaporized hydrogen peroxide steriliser: Amsco V-PRO1[®], Sakura Seiki) were measured by electrochemical detector (Polytron 7000[®], Dräger, sensitivity: 0-300ppm).

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
Method 2

- The HPV concentration on the plastic medical devices, experimental plastic panels (100 × 100 × 6mm), and the dummy of flexible fiberscope especially made were measured by the electrochemical detectors, sealing in the stainless steel containers.

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
Concentration of Hydrogen Peroxide Vapour in the Proximity of STERRAD®NX™ (A-D) and STERRAD® 200 (E)

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
Concentration of Hydrogen Peroxide Vapour in the Proximity of Amsco V-PRO1®

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About hydrogen peroxide sterilization

- Safety analyzer for hydrogen peroxide concentration in the proximity of each steriliser is not installed.
- And it is almost impossible to detect the concentration every time practically. So, any alarm apparatus as same as for ethylene oxide sterilisation is required for safety.

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
Hydrogen Peroxide Vapour from Plastic Test Panels after Low Temperature Hydrogen Peroxide Gas Plasma (LTHPGP) sterilisation (STERRAD[®] J&J)

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
Hydrogen Peroxide Vapour from PEI Panel (100D × 100W × 6 T mm)
Each Sealed in A Stainless Steel Container After Sterilization by STERRAD NX



Hydrogen Peroxide Vapour from Eleven Plastic Panels (100D × 100W × 6 T mm)
Each Sealed in A Stainless Steel Container After Sterilization by STERRAD NX®

- polyetherimide (PEI)
- polyethylene (PE)
- polytetrafluoroethylene (PTFE)
- polyamide 6 (PA6)
- polyamide 66 (PA66)
- polyethylene terephthalate (PET)
- polyetheretherketone (PEEK)
- thermoplastic polyurethane (TPU)
- polymethyl methacrylate (PMMA)
- polypropylene (PP)
- polycarbonate (PC)

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
Conclusion

- HPV sterilisation is the most useful procedure for the sterilisation of heat-labile medical devices at present. However, Safe and effective use of HPV steriliser should be re-evaluated.
- When HPV steriliser is employed, adequate air conditioning, donning gloves for bag handling, use of supplied-air respirator for door opening, adequate maintenance and other measures for the security of healthcare personnel should be seriously reevaluated.
- The CIs in the pouches with reusable items must be carefully checked before sterilization. It should be better to pack the devices with CIs just before sterilisation.
- For the safety of patients and hospital personnel, we must be careful for handling the medical devices sterilised by HPV and for the adverse events.

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